

iFDPacker simple User Guide


1. run iFDPacker.exe(by double click) in \Packer folder




2. Click Source [Browse...] button




Select the folder that Flashit.exe exist for pack DOS flash utility.




3. Click [Add] to select BIOS image



4. Select output file folder and specify a filename.(ex:mybios.exe).



5. Click [Create] button, will pop up a input arguments Dialog,
Input the argument(ex: /bios).
checked [Enable] to enable arguments filter for command line.
For example in below figure,
`C:>mybios.exe /nopause` (the /nopause argument will be passed into program)
`C:>mybios.exe /h` (the /h argument will not be passed into program)
Unchecked [Enable] all argument will be passed into program.
And click [OK] to close it.



6. If everything is OK, will pop up “Create package finished” Message Box, click [OK] to close it



Click [Exit] to close the application.

